


HL7 EASYCONNECT

A Service Based Interface Engine

Gain the benefits of an interface engine without the hassles!

It takes more than software to run an interface engine. Hospitals must develop the appropriate staff to implement and support critical interfaces 24 hours per day, 365 days per year. Maintaining your own interface engine can be an unreliable solution for IS Departments as employees turn-over, and illnesses and vacations threaten operational support.

At Iatric Systems, we have developed **HL7 EasyConnect**, which relieves your staff of having to maintain an interface engine. HL7 EasyConnect combines our robust interface engine software with our extensive implementation experience and our 24/7/365 support. With other vendors, you buy the software and YOU do all the work. With HL7 EasyConnect, you buy the software and WE do all the work. Major hospital systems with existing interface engines already have switched to HL7 EasyConnect in order to simplify their interface operations.



The screenshot shows the HL7 EasyConnect Remote Command Center interface. On the left is a tree view of interface configurations. The main area displays the 'Interface Router Status' table. The table has columns for No., Interface Name, Last Activity, Status, Queue, and Control. Row 13 is highlighted in yellow, indicating a disconnected interface.

No.	Interface Name	Last Activity	Status	Queue	Control
1	IN ADT from Meditech LIVE	1/7/2008 1:49:29 PM	Connected	---	Stop
2	IN ADT from Meditech TEST	1/3/2008 4:22:54 PM	Connected	---	Stop
3	OUT ADT to Witt	1/7/2008 1:49:31 PM	Connected	0	Stop
4	OUT ADT to Other	1/7/2008 1:49:31 PM	Connected	0	Stop
5	IN Reports from Witt	(No Activity Recorded)	Not Running	---	Start
6	IN Reports from Softmed	(No Activity Recorded)	Not Running	---	Start
7	OUT to Meditech NMI LIVE	(No Activity Recorded)	Not Running	0	Start
8	OUT to Meditech NMI TEST	(No Activity Recorded)	Not Running	0	Start
9	OUT ADT to Alpha TEST	1/3/2008 4:22:57 PM	Connected	0	Stop
10	OUT ADT to Alpha LIVE	1/7/2008 1:49:31 PM	Connected	0	Stop
11	OUT ADT to Intellidot TEST	1/3/2008 4:22:57 PM	Connected	0	Stop
12	OUT ADT to Intellidot LIVE	1/7/2008 1:50:01 PM	Connected	0	Stop
13	OUT ADT to OptiFlex TEST	1/7/2008 1:50:01 PM	Disconnected	69	Stop
14	OUT ADT to OptiFlex LIVE	1/7/2008 1:50:01 PM	Connected	0	Stop
15	OUT ADT/Orders to GE MUSE TEST	1/7/2008 10:28:14 AM	Connected	0	Stop
16	OUT ADT/Orders to GE MUSE LIVE	1/7/2008 1:49:29 PM	Connected	0	Stop
17	IN OE Orders from Meditech TEST	1/7/2008 10:28:11 AM	Connected	---	Stop
18	IN OE Orders from Meditech LIVE	1/7/2008 1:49:09 PM	Connected	---	Stop
19	IN Reports from GE MUSE TEST	1/7/2008 10:46:30 AM	Connected	---	Stop
20	IN Reports from GE MUSE LIVE	(No Activity Recorded)	Not Running	---	Start

HL7 EasyConnect enables data stream duplication and merging, field mapping, filtering, translation of data and more! For a demonstration or more information on **HL7 EasyConnect**, or any other Iatric Systems products or services, please contact us using the information below.