

One Device, Many Uses

Give your clinicians freedom and flexibility! Our **Mobile Solutions Suite** helps patient care providers accomplish all of the following activities from one easy-to-carry handheld device:

- Barcode specimen collection with **MobiLab**[®]
- Patient care documentation with **MobilCare**
- Barcode medication administration with **MobilMed**
- Charge entry with **MobilCharge**



Eliminate Mislabeled Specimens with MobiLab®

"We had no mislabeled specimens in any of the areas using MobiLab the very first month after go-live!"

Janet Johnson, BSN, MPH, RN-BC - Director of Nursing Informatics
Norman Regional Health System



MobiLab® can help improve patient safety at your hospital. By utilizing wireless handheld and printer technology, MobiLab is designed to eliminate patient identification and specimen labeling errors while increasing the productivity of your phlebotomy and nursing staff.

Positive Patient ID – Barcode technology is used to positively identify a patient by wristband barcode at the bedside. Users also are required to confirm a secondary identification method in order to complete the positive patient identification process.

Bedside Label Printing – Specimen collection labels are generated at the bedside through the use of portable label printers, significantly reducing the chance of mislabeled specimens.

Positive Specimen ID –

After specimen collection and labeling, barcode technology again is used to confirm that appropriate specimens have been collected and correct labels have been applied.

Motorola MC55



Zebra QL220 Plus Portable Printer



Key Advantages:

- Eliminates mislabeled specimens
- Improves stat turnaround times
- Decreases unnecessary venipunctures
- Runs on a variety of hardware – handhelds, laptops, PCs, mobile carts.

Ranked #1 in the KLAS

"Specimen Collection Barcoding" segment as published in the Top 20 Best in KLAS Awards report: 2007, 2008, 2009 and 2010.



(www.KLASresearch.com)

© 2010, KLAS Enterprises, LLC.

All rights reserved.)



Document Patient Care with MobilCare

Does your staff ever suggest that mobile carts are cumbersome to use? Now you can give them a convenient alternative! **MobilCare** gives clinicians the ability to view results and document patient care all from a handheld device.

MobilCare improves patient safety and ensures positive patient identification by requiring users to scan a patient's barcode wristband before documenting care. This true point-of-care solution enables staff to document care in a timely fashion at the patient's bedside.



Key Advantages:

- View all of your patients' interventions in next-due order
- View relevant patient data (lab results, allergies, demographics, I/Os, etcetera)
- Document interventions and assessments

Ensure the Five Rights of Medication Administration with MobilMed

MobilMed is a tool that can assist your staff in providing safe medication administration with the convenience of a handheld.

Positive Patient ID – Barcode technology is used to positively identify a patient by scanning the wristband barcode at the bedside. Users can be required to confirm a secondary identification method in order to complete the positive patient identification process.

Positive Medication ID – Barcode technology is used to positively identify the medication by scanning the NDC (National Drug Code) barcode number to ensure the right drug, the right dose, and the right route.

MobilMed enables clinicians to view eMARs (electronic Medication Administration Records) and document medication administration with ease. Scanning is built into the MobilMed handheld device and prevents the repeated scanner-to-keyboard switching necessary when using a PC or mobile cart.



Key Advantages:

- Ensures clinicians verify patients' "Five Rights"
- Enables customer defined screen (CDS) data entry during medication administration, including free text
- Can be used to replace PCs and mobile carts or as an adjunct to those devices



Verified patient and list of meds



Meds pending verification, dose instructions and label comments



Allergies displayed

Capture Patient Charges with MobilCharge

Streamline patient charge capture at the point of care. Using wireless handheld and barcode technology, **MobilCharge** enables staff to positively identify a patient via barcode wristband and scan chargeable items as they are used. MobilCharge is directly integrated with MEDITECH, so with the touch of a button a user can automatically update Materials Management (MM) and Billing/Accounts Receivable (B/AR) to reflect each patient item issued.

With MobilCharge, staff will capture previously lost charges and won't have to remember to key charges into the computer manually from sticky labels.



One Device, Many Uses – Mobile Solutions Suite

- Eliminate Mislabeled Specimens with **MobiLab®**
- Document Patient Care with **MobilCare**
- Ensure the Five Rights of Medication Administration with **MobilMed**
- Capture Patient Charges with **MobilCharge**

Iatric Systems' Mobile Solutions integrate directly with your MEDITECH system without the need for costly interfaces.

For more information or a demonstration, please contact us using the information below.



Ph: 978.805.4100 Fax: 978.805.4101 info@iatric.com iatric.com